



Date: - 15th December 2016

Corrigendum

for

Tender for Laparoscopic Instruments for the Department of General Surgery

NIT Issue Date	: 03 rd October, 2016
NIT No.	: Admn/Tender/59/2016-AIIMS.JDH
Pre-Bid Meeting	: 17 th October, 2016 at 03:45 PM.
Last Date of Submission	: 21 st November, 2016 at 03:00 PM.
Extended Last Date of Submission	: 03 rd January, 2017 at 03:00 PM.
Bid Opening	: 03 rd January, 2017 at 03:45 PM.

1. The following revised and additional specification will be added:-

1. Page No. 10, Annexure-I, Point No. 1:

For

Three part laparoscopic autoclavable KELLY Dissecting and Grasping Forcep, 360 degree rotational sheath, with connector pin for unipolar coagulation, size 5 mm, length 36 cm, long, double action jaws, with ergonomic plastic handle with larger contact area at the finger ring to avoid pressure sores, can be dismantled with the press of a button.

Read

Three part laparoscopic autoclavable KELLY/ **Maryland** Dissecting and Grasping Forcep, 360 degree rotational sheath, with connector pin for unipolar coagulation, size 5 mm, length 36 cm, long, double action jaws, with ergonomic plastic handle with larger contact area at the finger ring to avoid pressure sores, can be dismantled with the press of a button.

5 in number

2. Page No. 10, Annexure-I, Point No. 2:

For

Three part laparoscopic autoclavable CROCE-OLMI Grasping Forcep, 360 degree rotational sheath, size 5 mm, length 36 cm, long, double action jaws, ergonomic metal handle with plastic finger rings with larger contact area at the finger ring to avoid pressure sores, can be dismantled with the press of a button.

Read

Three part laparoscopic autoclavable CROCE-OLMI Grasping Forcep, 360 degree rotational sheath, size 5 mm, length 36 cm, long, double action jaws, ergonomic metal handle with plastic finger rings with larger contact area at the finger ring to avoid pressure sores, can be dismantled with the press of a button. 5 in number

3. Page No. 10, Annexure-I, Point No. 3:

For

Three part laparoscopic autoclavable Bowel Grasping Forcep, double action jaws, fenestrated, 360 degree rotational sheath, size 5 mm, length 36 cm, long, double action jaws, ergonomic metal handle with plastic finger rings with larger contact area at the finger ring to avoid pressure sores, can be dismantled with the press of a button.

Read

Three part laparoscopic autoclavable Bowel Grasping Forcep, double action jaws, fenestrated, 360 degree rotational sheath, size 5 mm, length 36 cm, long, double action jaws, ergonomic metal handle with plastic finger rings with larger contact area at the finger ring to avoid pressure sores, can be dismantled with the press of a button. . 5 in number

4. Page No. 10, Annexure-I, Point No. 4:

For

Three part laparoscopic autoclavable Right Angled Dissecting and Grasping Forcep, double action jaws, fenestrated, 360 degree rotational sheath, size 5 mm, length 36 cm, long, double action jaws, ergonomic metal handle with plastic finger rings with larger contact area at the finger ring to avoid pressure sores, can be dismantled with the press of a button.

Read

Three part laparoscopic autoclavable Right Angled Dissecting and Grasping Forcep, double action jaws, , 360 degree rotational sheath, size 5 mm, length 36 cm, long, double action jaws, ergonomic metal handle with plastic finger rings with larger contact area at the finger ring to avoid pressure sores, can be dismantled with the press of a button. .

5 in number

5. Page No. 10, Annexure-I, Point No. 5:

For

Three part laparoscopic autoclavable MANHES Grasping Forcep with ratchet, Tiger Jaws, 2X4 teeth, 360 degree rotational sheath, with connector pin for unipolar coagulation, size 5 mm, length 36 cm, long, double action jaws, with ergonomic plastic handle with larger contact area at the finger ring to avoid pressure sores, can be dismantled with the press of a button.

Read

Three part laparoscopic autoclavable MANHES Grasping Forcep with ratchet, Tiger Jaws, 2X4 teeth, 360 degree rotational sheath, with connector pin for unipolar coagulation, size 5 mm, length 36 cm, long, double action jaws, with ergonomic plastic handle with larger contact area at the finger ring to avoid pressure sores, can be dismantled with the press of a button.

5 in number

6. Page No. 10, Annexure-I, Point No. 6:

For

Three part laparoscopic autoclavable ANVIL Grasping Forcep, double action jaws, fenestrated, 360 degree rotational sheath, size 5 mm, length 36 cm, long, double action jaws, ergonomic metal handle with plastic finger rings with larger contact area at the finger ring to avoid pressure sores, can be dismantled with the press of a button.

Read

Three part laparoscopic autoclavable ANVIL/ Babcock Grasping Forcep, double action jaws, fenestrated, 360 degree rotational sheath, size 5 mm, length 36 cm, long, double action jaws, ergonomic metal handle with plastic finger rings with larger contact area at the finger ring to avoid pressure sores, can be dismantled with the press of a button. 5 in number.

7. Page No. 10, Annexure-I, Point No. 7:

For

Three part laparoscopic autoclavable MANHES Grasping Forcep with Hemostat style ratchet, duck-bill jaws, blunt, 360 degree rotational sheath, with connector pin for unipolar coagulation, size 5 mm, length 36 cm, long, single action jaws, with ergonomic plastic handle with larger contact area at the finger ring to avoid pressure sores, can be dismantled with the press of a button.

Read

Three part laparoscopic autoclavablev MANHES Grasping Forcep with Hemostat style ratchet, duck-bill jaws, blunt, 360 degree rotational sheath, with connector pin for unipolar coagulation, size 5 mm, length 36 cm, long, single action jaws, with ergonomic plastic handle with larger contact area at the finger ring to avoid pressure sores, can be dismantled with the press of a button. 2 in number

8. Page No. 10, Annexure-I, Point No. 8:

For

Three part laparoscopic autoclavable CLAW Grasping Forcep with ratchet, 2X3 teeth, 360 degree rotational sheath, size 10 mm, length 36 cm, long, single action jaws, with ergonomic metal handle, can be dismantled with the press of a button.

Read

Three part laparoscopic autoclavable CLAW Grasping Forcep with ratchet, 2X3 teeth, 360 degree rotational sheath, size 10 mm, length 36 cm, long, single action jaws, with ergonomic metal handle, can be dismantled with the press of a button. . 5 in number

9. Page No. 10, Annexure-I, Point No. 9:

For

Three part laparoscopic autoclavable Curved METZENBAUM Scissors, 360 degree rotational sheath, with connector pin for unipolar coagulation, size 5 mm, length 36 cm, long, double action jaws, with ergonomic plastic handle with larger contact area at the finger ring to avoid pressure sores, can be dismantled with the press of a button.

Read

Three part laparoscopic autoclavable Curved METZENBAUM Scissors, 360 degree rotational sheath, with connector pin for unipolar coagulation, size 5 mm, length 36 cm, long, double action jaws, with ergonomic plastic handle with larger contact area at the finger ring to avoid pressure sores, can be dismantled with the press of a button. 5 in number

10. Page No. 11, Annexure-I, Point No. 10:

For

Three part laparoscopic autoclavable Hook Scissors, 360 degree rotational sheath, with connector pin for unipolar coagulation, size 5 mm, length 36 cm, long, double action jaws, with ergonomic plastic handle with larger contact area at the finger ring to avoid pressure sores, can be dismantled with the press of a button.

Read

Three part laparoscopic autoclavable Hook Scissors, 360 degree rotational sheath, with connector pin for unipolar coagulation, size 5 mm, length 36 cm, long, double action jaws, with ergonomic plastic handle with larger contact area at the finger ring to avoid pressure sores, can be dismantled with the press of a button. 3 in number

11. Page No. 11, Annexure-I, Point No. 11:

For

Scissors Insert and outer tube, METZENBAUM Scissors, curved, length of blades 12 mm, double action jaws, size 5 mm, length 36 cm, package of 10, sterile, disposable.

Read

Scissors Insert and outer tube, METZENBAUM Scissors, curved, length of blades 12 mm, double action jaws, size 5 mm, length 36 cm, package of 10, sterile, disposable. 5 in number

12. Page No. 11, Annexure-I, Point No. 12:**For**

Scissors Insert with Outer Sheath, curved, double action jaws, spoon-shaped jaws, length of blades 17 mm, size 5 mm, length 36 cm, sterile, for single use, package of 10

Read

Read as: Scissors Insert with Outer Sheath, curved, double action jaws, spoon-shaped jaws, length of blades 17 mm, size 5 mm, and length 36 cm. 5 in number

13. Page No. 11, Annexure-I, Point No. 13:**For**

Surgical Sponge Holder, for atraumatic dissection of tissue layers, size 5 mm, length 30 cm, consisting of: Handle, Outer Sheath, insulated Sponge Holder Insert

Read

Surgical Sponge Holder, for atraumatic dissection of tissue layers, size 5 mm, length 30 cm, consisting of: Handle, Outer Sheath, insulated Sponge Holder Insert. 2 in number

14. Page No. 11, Annexure-I, Point No. 14:**For**

Ergonomic Plastic Handles with Larger contact area at the finger rings with connector pin for unipolar coagulation, with or without ratchet, compatible with above asked instruments.

Read

Ergonomic Plastic Handles with Larger contact area at the finger rings with connector pin for unipolar coagulation, with or without ratchet, compatible with above asked instruments
10 in number

15. Page No. 11, Annexure-I, Point No. 15:**For**

Dissecting-Electrode with exchangeable electrode L-tip, autoclavable, ergonomic handling, size 5 mm, length 36 cm, with connector pin for unipolar coagulation,

Read

Dissecting-Electrode with exchangeable electrode L-tip, autoclavable, ergonomic handling, size 5 mm, length 36 cm, with connector pin for unipolar coagulation. 3 in number

16. Page No. 11, Annexure-I, Point No. 16:**For**

Exchangeable Electrode Tip, L-Shaped, Autoclavable

Read

~~Deleted~~

17. Page No. 11, Annexure-I, Point No. 18:**For**

CUSCHIERI Liver Retractor, size 10 mm, length 36cm

Read

CUSCHIERI Liver Retractor, size 5 mm, length 36cm

18. Page No. 11, Annexure-I, Trocar & Cannula, Para No. 1:**For**

Cannula size: 13.5 **mm diameter**; should have multifunctional valve to prevent damage of sharp instruments and tip lens while passing through the cannula valve. It should have stopcock for CO₂ insufflation. Trocar should have pyramidal tip with pin holes near the tip for safety outlet of CO₂ gas. The working length of the cannula should be 115 mm.

Read

Cannula size: **13.5 mm diameter**; should have multifunctional valve to prevent damage of sharp instruments and tip lens while passing through the cannula valve. It should have stopcock for CO₂ insufflation. Trocar should have pyramidal tip with pin holes near the tip for safety outlet of CO₂ gas. The working length of the cannula should be 110 or more mm

19. Page No. 11, Annexure-I, Trocar & Cannula, Para No. 2:**For**

Cannula size: **11 mm diameter** with thread and rotating insufflations, should have multifunctional valve to prevent damage of sharp instruments and tip lens while passing through the cannula valve. It should have stopcock for CO₂ insufflation. The working length of the cannula should be 105 mm.

Read

Deleted

20. Page No. 11, Annexure-I, Trocar & Cannula, Para No. 3:**For**

Telescope Stopper, size 11 mm

Read

Deleted

21. Page No. 12, Annexure-I, Trocar & Cannula, Para No. 4:**For**

Cannula size: **6 mm diameter** with thread and rotating insufflations, should have multifunctional valve to prevent damage of sharp instruments and tip lens while passing through the cannula valve. It should have stopcock for CO₂ insufflation. The working length of the cannula should be 105 mm.

Read

Deleted

22. Page No. 12, Annexure-I, Trocar & Cannula, Para No. 5:**For**

Telescope Stopper, size 6 mm

Read

Deleted

23. Page No. 12, Annexure-I, Trocar & Cannula, Para No. 7:**For**

Cannula size: **6mm diameter**; should have multifunctional valve to prevent damage of sharp instruments and tip lens while passing through the cannula valve. It should have stopcock for CO₂ insufflation. Trocar should have pyramidal tip with pin holes near the tip for Safety outlet of CO₂ gas. The working length of the cannula should be 105 mm.

Read

Cannula size: **5.5-6mm diameter**; should have multifunctional valve to prevent damage of sharp instruments and tip lens while passing through the cannula valve. It should have stopcock for CO₂ insufflation. Trocar should have pyramidal tip with pin holes near the tip for Safety outlet of CO₂ gas. The working length of the cannula should be 105 mm.

24. Page No. 12, Annexure-I, Trocar & Cannula, Para No. 8:**For**

Double reducer, 13.5/10 & 13.5/5mm

Read

Double reducer, 13.5/10 & 13.5/5mm. 10 number each

25. Page No. 12, Annexure-I, Trocar & Cannula, Para No. 9:

For

Reducer, 11/5 mm

Read

10 number each

26. Page No. 12, Annexure-I, Veress Needle, Para No. 1:

For

Veress Pneumoperitoneum needle with spring loaded blunt inner cannula, LUER-Lock, autoclavable, diameter 2.1mm , working length 10- 18cms.

Read

VeressPneumoperitoneum needle with spring loaded blunt inner cannula, LUER-Lock, autoclavable, diameter 2.1mm, working length 10- 18cms. 10 in number

27. Page No. 12, Annexure-I, Injection Needle, Para No. 1:

For

Injection needle, LUER-lock, diameter 1.2 mm, size 5 mm, length 36 cm.

Read

Injection needle, LUER-lock, diameter 1.2 mm, size 5 mm, length 36 cm. 3 in number

28. Page No. 12, Annexure-I, after Injection Needle:

For

Suction and irrigation 2 way system.

Read

Simple 2 way up and down suction cannula with 36 cm. 5 mm. 5 in number

29. Page No. 12, Annexure-I, after Injection Needle: LAPAROSCOPIC NEEDLE HOLDERS, Para No. 1:

For

Macro needle holder with tungsten carbide insert, dismantling into three parts namely, outer tube, handle and inserts, ergonomic pistol handle, with disengageable ratchet, jaw curved to left, size: 5mm, length: 33cm for use with suture material size: 0/0 to 7/0.

Read

Macro needle holder with tungsten carbide insert, dismantling into three parts namely, outer tube, handle and inserts, ergonomic pistol handle, with disengageable ratchet, jaw curved to left, size: 5mm, length: 33cm for use with suture material size: 0/0 to 7/0.

Right curved tip.3 in number

Left curved tip. 2 in number

One coaxial. 1 in number

30. Page No. 12, Annexure-I, after Injection Needle: LAPAROSCOPIC NEEDLE HOLDERS, Para No. 2:

For

Macro needle holder with tungsten carbide insert, dismantling into three parts namely, outer tube, handle and inserts, ergonomic pistol handle, with disengageable ratchet, jaw curved to right, size: 5mm, length : 33cm for use with suture material size: 0/0 to 7/0.

Read

Deleted

31. Page No. 12, Annexure-I, after Injection Needle: Laparoscopic Bipolar Forceps, Para No. 1:

For

Three part laparoscopic Bipolar Coagulating Grasping Forceps with overload tissue protection, CLERMONT- FERRAND MODEL, Wide jaws, with connector pin for bipolar coagulation, 360 degree rotational sheath, size 5 mm, length 36 cm, long, single action jaws,

ergonomic Plastic handle with larger contact area, O can be dismantled with the press of a button.

Read

Three part laparoscopic Bipolar Coagulating Grasping Forceps with overload tissue protection, Wide jaws, with connector pin for bipolar coagulation, 360 degree rotational sheath, size 5 mm, length 36 cm, long, single action jaws, ergonomic Plastic handle with larger contact area, O can be dismantled with the press of a button. 2 in number

32. Page No. 12, Annexure-I, after Injection Needle: Laparoscopic Bipolar Forceps, Para No. 2:

For

Bipolar HF Connecting Cable

Read

Bipolar HF Connecting Cable: 3 in number

33. Page No. 12, Annexure-I, after Injection Needle: Laparoscopic Bipolar Forceps, Para No. 3:

For

Clip Applicator: 01

Medium Large clip applicator dismantable rotating size: 10mm, length 36cm, for Titanium clips with ratchet to lock the jaw holding the clip.

Read

Medium Large clip applicator dismantable rotating size: 10mm, length 36cm, for Titanium clips with ratchet to lock the jaw holding the clip. 2 in number

34. Page No. 12, Annexure-I, after Injection Needle: Laparoscopic Bipolar Forceps, Para No. 4:

For

Clip Applicator: 01

Large clip applicator dismantable rotating size: 10mm, length 36cm, for Titanium Clips with ratchet to lock the jaw holding the clip.

Read

Large clip applicator dismantable rotating size: 10mm, length 36cm, for Titanium Clips with ratchet to lock the jaw holding the clip. 3 in number

35. Page No. 12, Annexure-I, after Injection Needle: Laparoscopic Bipolar Forceps, Para No. 4:

For

Titanium Clips: 10

Titanium clips medium large & large, box with 16 sterile cartridges, 10 clips each For use with clip applicator.

Must supply suitable Autoclavable plastic double tray for sterilization and storage of at least 20 hand instruments preferably from OEM 1 in number.

36. Page No. 12, Annexure-I, after Injection Needle: Laparoscopic Bipolar Forceps, after Para No. 3:

Add items with specification

10 mm Hassan's cannula. Used for TEP repair. 1 in number

All instrument should have European CE and USFDA approved

**Administrative Officer
AIIMS, Jodhpur**