

All India Institute of Medical Sciences, Jodhpur

No. Admn/Prop/63/2018-AIIMS.JDH

Dated 31th January, 2019

Subject: Proposal for Purchase of reagents for Haematology Analyzer System Model XN 1000B₃ on Proprietary basis at AIIMS Jodhpur – Inviting Comments/Objection thereon.

The Institute is in the process of purchasing reagents for Haematology Analyzer System Model XN 1000B3 on Proprietary basis at AIIMS, Jodhpur from M/s Transasia Bio Medicals Ltd., Mumbai through authorized distributors M/s Suncity Medicare, Jodhpur. The proprietary certificate submitted by the user department are attached & uploaded on website.

The above documents are being uploaded for open information to submit objection, comments, if any, from any manufacturer regarding proprietary nature of the above said solutions within 7 days from the date of issue/uploading of the notification. The comments should be sent to administrative officer, AIIMS Jodhpur on or before 06th February, 2019 up to 03:00 PM, failing which it will be presumed that any other vendor is having no comment to offer and case will be decided on merits.

Administrative Officer

Enclosure: Related documents enclosed.



Transasia Bio-Medicals Ltd., Transasia House, 8 Chandivali Studio Road, Andheri (East), Mumbai - 400 072 Tel: +91 22 4030 9000 Fax: +91 22 2857 3030 Email : transasia@transasia.co.in CIN : U33110MH1985PLC036198



Monday, 28 January 2019 Ref. TBM/AIIMS/18-19/902

To, Administrative Officer All India Institute of Medical Sciences, Basni, Ilnd Phase, Jodhpur

Subject: Proprietary Certificate for Sysmex 6 part differential Hematology Analyzer **Model**-XN 1000 B3 at AIIMS, Jodhpur

Dear Sir,

This is to certify that the following reagents have to be used on our Sysmex 6 part differential Hematology Analyzer Model-XN1000 B3. The reagents are imported from Sysmex Corporation. This is a closed system Instrument and no other Reagents can be used on this Instrument except the following:-

- CELL PACK DCL (20L)
- SULFOLYSER (1.5L x 2)
- Lysercell (5Ltr) (WNR)
- Fluorocell (82mL* x 2) (WNR)
- Lysercell (5Ltr) (WDF)
- Fluorocell (42mL* x 2) (WDF)
- Fluorocell Plt (12mL* x 2)
- CELLPACK DFL (1.5L x 2) (RET)
- Fluorocell (12mL* x 2) (RET)
- Cell Clean (50ml)
- QC XN Check L/H/N (3x4.5ml)

We cannot guarantee the accuracy and precision of our analyzers if any other reagents are used on Sysmex XN1000 B3.

Thanking You

FOR TRANSASIA BIO MEDICALS LTD.

Ram J Kajla 28-1.2019

(Regional Sale Manager) Ph. 9309440143 Date: 28-01-2019

www.transasia.co.in






Transasia Bio-Medicals Ltd., Transasia House, 8 Chandivali Studio Road, Andheri (East), Mumbai - 400 072 Tel: +91 22 4030 9000 Fax: +91 22 2857 3030 Email : transasia@transasia.co.in CIN : U33110MH1985PLC036198



Monday, 28 January 2019 Ref. TBM/AIIMS/18-19/901

To, Administrative Officer All India Institute of Medical Sciences, Basni IInd Phase, Jodhpur



Subject: RE: Sysmex Reagent Quotation No. HLL/PCD/PMSSY/AIIMS-II/19/15-16/4701 dt. 29-Dec-15-6 nos. XN 1000 B3 from AIIMS, Jodhpur

Dear Sir,

This has reference to NOA No. HLL/PCD/PMSSY/AIIMS-II/18/15-16/4701 dt. 29-Dec-15 for XN 1000 B3 supplied to AIIMS, Jodhpur Below are the Rates of Proprietary Sysmex Reagents and Consumables for XN1000 B3:-

Product Name CELL PACK DCL (120967)		Pack Size
WNR	Lysercell (121531)	5 Ltr
	Fluorocell (120975)	82mL* x 2
WDF	Lysercell (121530)	5 Ltr
	Fluorocell (120977)	42mL* x 2
RET	CELL PACK DFL (120980)	1.5L x 2
	Fluorocell (120981)	12mL* x 2
	Fluorocell PLT (120982)	12mL* x 2
	Cell clean (120515)	50 ml
	QC – XN Check L/H/N	3x4.5 ml

*GST Extra

Above Sysmex Make Reagents & Consumables will be supplied through our Authorized Distributor M/s Suncity Medicare, 4 Happy Home Tower, 4th C Road, Sardarpura, Jodhpur-342001 (Rajasthan).

Note: This Rate Quotation is Valid upto 31/03/2019 and GST Extra on above Price. Supply period will be 30 to 45 Days from Confirm Purchase Order Date.

Thanking You

FOR TRANSASIA BIO MEDICALS LTD.

Ram J Kajla 28.1.2019

(Regional Sale Manager) Ph. 9309440143

-Jay Six

Date: 28-01-2019

