



Date: - 11<sup>th</sup> January, 2019

Corrigendum  
For  
Office Hysteroscopy Set for the Department of  
Obstetrics and Gynaecology

NIT Issue Date	: 12 <sup>th</sup> October, 2018
NIT No.	: Admn/Tender/39-2/2018-AIIMS.JDH
Pre-Bid Meeting	: 22 <sup>nd</sup> October, 2018 at 03:15 PM
Earlier Last Date of Submission	: 14 <sup>th</sup> January, 2019 at 03:00 PM
Extended Last Date of Submission	: 30 <sup>th</sup> January, 2019 at 03:00 PM
Bid opening	: 31 <sup>st</sup> January, 2019 at 03:15 P.M

The following revised and additional specification will be added:-

1. Page No. 11, Point 1, bullet no. 1

For

Office Hysteroscopy 30° (Diameter 1.9-2.7mm)- 1 Nos

Read

Office Hysteroscopy 30° (Diameter 1.9-2.9 mm)- 1 Nos

2. Page No. 11, Point 1, bullet no. 2

For

Diagnostic sheath must not exceed 3.6 mm when assembled for single flow.

Read

Diagnostic sheath must not exceed 3.6 mm when assembled for single flow. 1 Nos

3. Page No. 11, Point 1, bullet 3:

For

Therapeutic sheath must not exceed 4.5mm when assembled for continuous flow with 5fr working channel- 1 Nos

Read

Therapeutic sheath must not exceed 4.5mm when assembled for continuous flow with 5fr working channel- 1 Nos

4. Page No. 11, "Camera Head", Point 2, bullet 3:

For

Signal noise: 50 dB, 1 Nos

Read

Signal to noise ratio: 50 dB, 1 Nos

**5. Page No. 11, Point 3, “Monitor” bullet 2:**

**For**

It should have USB port for Capturing videos and Stills in External USB drive

**Read**

It should have USB port for Capturing videos and Stills in External USB drive